
# STRAIGHT & TILTED

CE  
0123

SURGICAL GUIDE **REF** 4880

IT

EN

DE

FR

ES

PT

DA

NL

SV

EL

ZH

HR

NO

LV

LT

PL

RO

HU

CS

TR

BG

RU

SL

KO

KK

UA

Rev. 01 – January 2020


META

 **C.G.M. S.p.A.**  
**DIVISIONE MEDICALE META**

Via E. Villa, 7 - 42124 Reggio Emilia - ITALY  
Tel.: +39 0522.50.23.11 - Fax: +39 0522.50.23.33  
[www.metahosp.com](http://www.metahosp.com)

**FIG. 1**

**FIG. 2**

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- The device is supplied non-sterile. The device must be cleaned, disinfected and sterilized before each use.
- The device must be used exclusively by skilled medical staff.
- The surgeon must assess the suitability of patients for treatment and the appropriate surgical procedure.
- The surgeon must take all necessary precautions to prevent the risk of infection and complications associated with the surgical procedure.
- Do not use the device if it is worn or damaged
- Do not use thermal disinfection methods.
- Do not use automatic disinfection methods.
- Clean, disinfect and sterilize the device before use.
- Always wear sterile gloves when handling the device and strictly follow the procedures in order to guarantee a sterile environment.
- Discard after use in special medical waste containers in compliance with the regulations in force.
- Meta does not shoulder any r

### 1 – INTENDED USE

The "Straight & Tilted Surgical Guide" provides the surgeon with a reference regarding distance and angle of the hole to be prepared in the "ALL ON FOUR" technique in oral surgery.

The guide is supplied non-sterile and tested for 30 uses in conditions of normal use and maintenance and in accordance with these instructions for use.

It is recommended to keep a record of the uses so as not to exceed the recommended number of uses.

The "Straight & Tilted Surgical Guide" has the following advantages: prevents reflections caused by any lamps thanks to the fact that it is made of 316L stainless steel with matt satin finish; provides the surgeon with an ideal reference point for the pillar thanks to the pre-angle achieved by marking the inclinations (0°, 17° and 30°); allows blocking the angle of the guide after modeling on the patient thanks to an adjustment screw.

In any event, the guide is only intended as a reference. The surgeon must assess case by case the best possible angle and distances of the pillars. Useful tips:

In the case of the lower arch, it is important to identify the mental foramen and inferior alveolar nerve. The final position of the implant should be before the foramen avoiding the loop of the nerve.

In the case of the upper arch, the position of the implant should avoid the maxillary sinus.

### 2 – BEFORE AND AFTER USE

Before each use:

always examine the device in order to check for any wear or damage; perform the cleaning, disinfection and sterilization procedures described in steps 3, 4 and 5.

Immediately after the surgical procedure:

perform the cleaning, disinfection and sterilization procedures described in steps 3, 4 and 5.

### 3 – CLEANING THE DEVICE

Cleaning consists of the following steps:

(a) The instrument must be subjected to a pre-cleaning phase as follows:

- disassemble the instrument as shown in Figures 1 and 2;
  - place the devices in a container separately from other materials to prevent damage or contamination;
  - clean the components with deionized cold water (15°C – 18°C) using a soft-bristle brush to remove any residue;
  - dry well using a paper or cloth towel with low particle release.
- Check for any damage.

After manual cleaning proceed with ultrasonic cleaning.

(b) Ultrasonic cleaning should always be carried out and is a manual but not a replacement treatment.

For ultrasonic cleaning:

- place the device in a container separately from other materials to prevent damage or contamination;

- it is recommended to use a detergent with the following characteristics: non-ionic surfactant;

- at least based on quaternary ammonium salts;

- with proven bactericidal, yeasticidal and mycobactericidal efficacy; suitable for use in an ultrasonic cleaner.

The manufacturer has validated the following products and methods:

Gigasept Instru AF. The solution was prepared, following the manufacturer's instructions, at a concentration of 2.5% of the active ingredient belonging to the category "quaternary ammonium salts", with water at a temperature of 30°C, contact time 5 minutes, in an ultrasonic cleaner at 40 KHz and power 140 W;

Surfanios. The solution was prepared, following the manufacturer's instructions, at a concentration of 0.0625% of the active ingredient belonging to the category "quaternary ammonium salts", with warm water (30°C - 32°C), contact time 5 minutes, in an ultrasonic cleaner at 40 KHz and power 140 W.

Other detergents and methods must be validated by the user.

For the methods of use (concentration of the solution, temperature and contact time) always refer to the manufacturer's instructions of each detergent to obtain the required microbicidal efficacy for reprocessing.

The use of a specific ultrasonic cleaner for dental surgeries is recommended. Make sure that the surface of the device components is completely covered, including the hollow parts, in order to allow the solution to produce its effects. Do not use oxidizing acids such as sulfuric acid, nitric acid.

After ultrasonic cleaning, rinse thoroughly for at least 3 minutes using sterile distilled water or deionized water to completely remove the residues.

The cleaning solution must be replaced after each use.

For the cleaning process it is recommended to follow the requirements of ISO 15883-1.

Dry using a paper or cloth towel with low particle release before proceeding with the subsequent disinfection phase.

### 4-DISINFECTION

For chemical disinfection fully immerse the disassembled device in sterile distilled water or deionized water and disinfectant.

It is recommended to use a disinfectant with the following characteristics:

- at least based on quaternary ammonium salts with a concentration of 5%;
- with proven bactericidal, yeasticidal and mycobactericidal efficacy.

The manufacturer has validated the following product and method:

I-SAN HYGIENE SICURA: The solution was prepared undiluted, contact time 15 minutes, following the manufacturer's instructions for bactericidal efficacy. The concentration of active ingredient belonging to the category "quaternary ammonium salts" in this solution is 5%.

Other disinfectants and methods must be validated by the user.

For the methods of use (concentration of the solution, temperature and contact time) always refer to the manufacturer's instructions of each detergent to obtain the required microbicidal efficacy for reprocessing.

Make sure that the surface of the device components is completely covered, including the hollow parts, in order to allow the solution to produce its effects. Do not use oxidizing acids such as sulfuric acid, nitric acid. Do not use thermal disinfection methods. Do not use automatic disinfection methods because they could damage the device. The device has only been tested for manual disinfection. Meta is not liable for damage to devices disinfected using invalid methods.

After disinfection rinse thoroughly for at least 3 minutes, using sterile distilled water or deionized water to completely remove the residues. The disinfectant solution must be replaced after each use.

Proceed to the drying phase through the following phases:

- use a paper or cloth towel with low particle release;
- dry in an oven at a temperature of 35°C for at least 1 hour.

The drying phase is particularly important as the introduction of not perfectly dry material into autoclaves can compromise the subsequent sterilization process.

The cleaning, washing, disinfection and drying procedures before and after using the medical device are intended to remove post-surgery organic

residues and lower the bioburden of the device. The surgeon should periodically ensure that the bioburden of the device does not exceed the alert limit for the effectiveness of the subsequent sterilization process. After this phase, go ahead to sterilization.

## 5 – STERILIZATION OF THE DEVICE

All the components of the device can be sterilized in steam autoclaves. During sterilization, protect the disassembled device from contact with other substances or materials in order to avoid any alterations, sealing it in a special bag - consisting of medical paper and polyester/polypropylene film - and proceed with sealing it as indicated by the bag manufacturer.. The thermal sealing process must be validated.

The polymer material of the bag must comply with DIN EN ISO 11607 Parts 1 and 2.

The paper must comply with DIN EN ISO 11607 Part 1 and DIN EN 868 Part 3.

The bag must have a Class 1 chemical colour change indicator on the outside in compliance with UNI EN ISO 11140-1, indicating whether the package has been sterilized.

Sterilization must be performed in class "B" steam autoclaves in accordance with UNI EN 13060, which assures the sterility of hollow and/or porous bodies.

The autoclave must undergo regular maintenance as indicated by the manufacturer in the user manual.

The steam autoclave (wet heat) sterilization method used by META complies with the ISO 17665 series standards. It is recommended to sterilize the device with the following parameters:

- temperature  $121^{\circ}\text{C} \pm 3^{\circ}\text{C}$ ;
- time 20 minutes;
- pressure 1 bar  $\pm 0.5$  bar;
- presence of saturated steam.

The sterilization process validated by the manufacturer guarantees a sterility level (SAL) of  $10^{-6}$ .

In any case, refer to the autoclave manufacturer's indications based on the specific load configuration of each dental surgery and the characteristics of the device to be processed.

It is recommended to always use sterilization cycles validated in accordance with ISO 17665-1 and to respect the validated loading configuration.

For autoclave loading operations, it is important to follow a few simple rules: do not overlay the material by stacking it;

do not obstruct the exhaust of the autoclave;  
do not cover the area of the filters and valves of the baskets;  
prevent the bags from touching the inner walls of the sterilizing chamber;  
the bags should be arranged by "size" (in a position parallel to the steam flow);

lay the bags one next to the other without pressing them and with the porous part facing upwards, using the dedicated perforated trays to contain them;  
place hollow and heavy objects on the lower shelf;  
do not overload the autoclave.

During the process it is necessary to record the temperature, the pressure of the saturated steam and the times in a dedicated logbook.

### 5.1 - TESTS TO BE PERFORMED

We recommend you perform the following tests on the autoclave for steam sterilization.

**VACUUM TEST:** daily test at the start of the sterilization before starting the sterilization cycles.

**HELIX TEST:** daily test after the vacuum test and before normal use of the autoclave.

With each load you run, it is recommended to insert a **CHEMICAL STERILITY INDICATOR** inside a package, in the same conditions as the material being sterilized. This test measures the parameters: temperature, pressure and exposure times. In addition, the bags themselves must have a **CHEMICAL PROCESS INDICATOR** in order to check whether or not the package has been exposed to a sterilization cycle.

**BIOLOGICAL INDICATOR:** it is recommended to carry out this test at least once a year.

### 5.2 - UNLOADING THE STERILE MATERIAL

At the end of the sterilization cycle you have to be careful that there is no moisture residue on the material itself, therefore, place it on a dry and clean

surface protected from air movements, covering it with a clean and dry cloth for at least 10 minutes.

### 5.3 - STORAGE OF THE STERILE MEDICAL DEVICE

Once treated, the medical device can be stored for up to 6 months in a drawer or a closed place to prevent any actions that may damage it, such as bumps, dropping the packs, rubbing. Indicate the sterilization date and its expiry on the bag. Upon expiry, steps 3, 4 and 5 must be repeated. Check the integrity of the bag before use.

### 5.4 - TRACEABILITY OF THE STERILIZATION PROCEDURE

The traceability system is a recommended permanent recording system that through a series of data indicated on the package allows you to identify and possibly track the sterilized device at any time of its journey (from unloading it from the autoclave to its use) and to trace back to the operator responsible for the process and the patient. It is recommended to indicate the following data on each package:

- acronym of the operator responsible for the sterilization cycle;
- number of the sterilization cycle;
- packing and expiration date.

### 5.5 - LOGBOOK

It is recommended that the recordings be stored in a specific document such as "sterilization log" or "sterilization book". The log should show: the result of the vacuum test (printing or recording in software) or hand-transcription of the temperature, pressure and time parameters; the Helix test;

printout of the process parameters (temperature, pressure, time) of each cycle, otherwise handwritten recording; the date and identification code of the cycle.

All the sterilization process documents must be archived and available in the workplace.

### 6 – DEVICE RE-USE

Before any re-use, check that the sterile bag is intact, that sterilization has not expired, that the device is clean, free of dust, dirt and any bone residue. If not, repeat the cleaning, disinfection and sterilization phases as described in steps 3, 4 and 5. Also check that the device is not damaged, in which case do not use it and discard it in the appropriate medical waste containers.

These instructions for use as well as further information about the instrument are available on the website: [www.metahosp.com](http://www.metahosp.com)

Please take the time to read this information to prevent any errors and ensure good operation of the product.

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