

Declaration of Conformity

GKE Seal-Test

**for daily routine control of sealing seams
and for the operational qualification of different sealing machines
in accordance with EN ISO 11607-2**

Art. No.*	Product code	Quantity	Dimension	Application	Standard
200-311	SEAL-TEST	100	173 x 76.2 mm	Seal indicator to check paper/plastic film-pouches in rotary sealing machines	EN 868-3
200-312		250			
200-321	HDPE-SEAL-TEST	100	173 x 76.2 mm	Seal indicator to check Tyvek® pouches in rotary sealing machines	EN 868-9
200-331	BAR-SEAL-TEST	100	250 x 76.2 mm	Seal indicator to check paper/plastic film pouches in bar sealing machines	EN 868-3

(*) On all GKE packages, an additional letter code has been added to the 6-digit article number. The additional letter code refers to the language and/or customized version. It is only added on the outside label, the inside of the pack is identical to the article numbers and the above table. All articles with the same 6-digit number have the same specifications.

The above-mentioned seal-tests are manufactured according to the corresponding standards and therefore conform in general with the standard requirements, which are not changing from batch to batch. In contrast to biological indicators with batch-related modifications, an individual batch-related certificate does not make sense because the specifications of the seal-tests do not change batch by batch.

We hereby declare under our sole responsibility, that the above-mentioned seal-tests have been tested in the GKE application laboratory and are in conformity with the above-mentioned European Standard.

This document certifies that the above performance criteria and the GKE test requirements for quality control are met. The continuous quality is guaranteed by our quality management system according to EN ISO 13485*.

Waldems, 2022-04-27



Dipl.-Ing. Dr. Ulrich Kaiser
R & D-Manager

* This certificate is available on the GKE homepage www.gke.eu.